


fp fp fp fp fp fp fp fp fp fp fp fp fp fp fp fp

Instructions Floor Pressure Mat - PMF

The Floor Pressure Mat alerts to foot pressure.


Moments after someone steps onto the mat the alarm will sound. Mats are often positioned in a doorway, hallway, or next to a bed.



PLUG-MATCHED or PAGER-LINKED
Supplied with a 3m lead ready to plug into your existing call system or our pager system.



ANTI-SLIP INCLUDED
-Grips to any surface
-Keep mat firmly in position
-Machine washable



OPTIONAL BARRIER MAT
-Protects mat from wear
-Significantly prolongs mat life
-Thick heavy duty rubber
-Absorbant cotton carpeted layer
-Keep mat firmly in position
-Machine washable





PAGER-LINKED

If you have a **Pager-Linked** floor pressure mat you will have a transmitter box included. Plug the mat into the socket on the box.



PLUG-MATCHED

Plug your floor pressure mat into your existing call system using the matched plug for the connection on your call point. An adapter lead may be supplied depending on the call system type.



SPLITTER PLUG

An optional splitter plug allows an existing call bell to be used in the same socket as your new floor pressure mat. The splitter converts a socket in a call system so it can take two inputs.

Mat Care

When the mat is not in use for long periods of time it is best not to fold the mat.

Never position a chair leg or sharp/heavy object on the mat as this can damage its internal workings.

The mat can be cleaned with a damp cloth using normal cleaning products or alcohol.

Safety

Remember, daily system tests should be carried out to ensure correct functioning of the unit. Usage should be incorporated within safety manuals and procedures. Range tests should be carried out at least once a week, more often if critical criteria apply. This should involve testing the unit past its required range. If the unit has been dropped or it is worn by a person involved in an accident the unit should be tested again before re-use.

Cleaning

The unit can be cleaned using a damp cloth and a small amount of disinfectant, alcohol based cleaning product or diluted chlorine-based cleaning solution. Do not submerge as the unit is not waterproof.

Care

DO NOT subject this equipment to: Mechanical shock, Excessive humidity, Extremes of temperatures, Corrosive Liquids.

This equipment is designed primarily for indoor use and is not water resistant. It must not be used in classified hazardous areas including areas containing explosive or flammable vapours. Consult your local product dealer for further information.

Specifications:

Dimensions(mm)	560mm X 720mm X 3mm
Weight	400g

Liability

Frequency Precision does not accept any liability for any damage or injury, howsoever caused as a result of misuse of this equipment. It is the responsibility of the user to ensure that the equipment is operated in the manner for which it was intended and that it is the correct item of equipment for the required task.

All systems can fail and it is the responsibility of the user to carry out regular tests and to determine the suitability of this equipment for any application.

Repair and replacement

Frequency Precision will refund payment for any unit returned within 30 days of purchase as unsuitable for the intended purpose. Un-damaged units will be repaired free of charge within the first 12 months.

Literature

Frequency Precision Ltd operates a policy of continual improvement and therefore reserves the right to modify and change any specification without prior notice.

While every possible care has been taken in the preparation of its manual, we do not accept any liability for the technical or typographical errors or omissions contained herein, nor for incidental or consequential damages arising from the use of the material.

Disposal

At the end of the working life of the product it must not be disposed of with household waste but returned to Frequency Precision Ltd or disposed of at a collection point for the re-cycling of electrical and electronic equipment.