

John Preston & Co



John Preston

*Bariatric Equipment
Package*



*Provision & Dignity should
come in all shapes and
sizes!*

Contents

Specialist Seating

Page 3 – Shannon High back chair

Page 4 – Teal 747 High Back chair

Page 5 – Dalerise Riser Recliner

Page 6 – Preston Rimini Riser Recliner

Walking Aids

Page 7 – Four Wheeled Rollator

Page 8 – Walking Frame

Page 9 – Folding Walking Frame

Page 10 – Wheeled Walking Frame

Moving & Handling Aids

Page 11 – Transfer Board

Page 12 – Rotunda turning Aid

Page 13 – Lina Mobile Hoist

Page 14 – AirSlide Lateral Transfer system

Page 15 – Porta 200 Portering chair

Bathroom Aids

Page 16 – Titan Shower Commode

Page 17 - Titan Plus Shower Commode

Page 18 – Zenith Static Commode

Page 20 – Terms and Conditions

Specialist Seating

Shannon High Back Chair



MAX USER WEIGHT 40ST/255KGS

Seat width – 24”

Seat depth – 19.5”

Seat height adj – 19”-23”

The Shannon Bariatric Chair is an excellent choice for larger mobile users who require basic seating and support.

Teal 747



MAX USER WEIGHT 40ST/255KGS

Seat width – 24”

Seat depth – 18.5”

Seat height adj – 17”-20”

This bariatric high back chair is an excellent choice for larger mobile users who require basic seating and support. Comes with a built in high risk pressure relief cushion as standard. Optional features include drop side arms, arm pads and height adjusters.

Dalerise Riser Recliner



MAX USER WEIGHT 28st

Seat width – 20”

Seat depth – 18”

Seat height – 18”

This Riser Recliner is an excellent choice for larger mobile users who require basic seating and support and a range of positions. A heavy duty triple motor system with easy to use hand-control; as well as a very stable anti-entrapment frame and battery back-up facility as standard. This chair is covered in vinyl and is easily cleaned for infection control purposes.

Preston Rimini



MAX USER WEIGHT 40st

Seat width – 26”

Seat depth – 18-22”

Seat height – 18”

This Riser Recliner is an excellent choice for very large users who require long term seating, high levels of all over pressure relief and good support.

This chair offers users a wide range of positions with a superb level of comfort. A heavy duty dual motor system with easy to use hand-control and battery back-up facility as standard. This Preston Rimini is covered in multi stretch breathable and wipeable fabric on patient contact surfaces and is easily cleaned for infection control purposes. Wooden hand knuckles provide improved grip point for transfers. This chair is excellent for users with heavy limbs and the fully adjustable waterfall backrest allows perfect contouring.

Walking Aids

Four - Wheeled Rollator



MAX USER WEIGHT 50ST/315kgs

Seat width – 18.5”

Seat depth – 14”

Seat height – 21.5”

High seat to floor allows easier transfer for user. Extra wide padded seat and backrest. Adjustable height push handles, folds for storage and transportation. Locking brakes. Puncture proof tyres. Storage basket.

Heavy Duty Walking Frame



Max User weight 318kgs/ 40ST

Heavy Duty Steel Construction

Comfortable moulded hand grips

Non-slip rubber feet

Height adjustable

Maximum width – 30”

Bariatric Folding Walking Frame



Max User weight 200kgs/ 40ST

Comfortable moulded hand grips

Non-slip rubber feet

Easily folded for transportation

Bariatric Wheeled Walking Frame

Height adjustable 78-101cm

Maximum width – 30"



Max User weight 325kgs/ 42ST

Comfortable moulded hand grips

Non-slip rubber feet

Easily folded for transportation

Height adjustable 78-101cm

Maximum width – 30”

Moving & Handling Aids

Transfer Board



Max User weight – 285kgs

The Q2 is a shatter proof polypropylene board impregnated with antimicrobial agent for improved infection control.

Smooth finish for ease of transfer.

Rotunda



Max User weight 70st

Robust Heavy duty standing

Turns and locks in 4 positions

Rubber handgrips

Rear wheels for moving frame when not in use

Lina Mobile Hoist



Max User Weight 200kgs

Mobile hoist with 2 point spreader bar

Compact design

Powered lift and powered opening legs

Easy to use handset

Emergency lowering facility

Locking castors

Air Slide



Max User weight 544kgs

The Airslide is inserted under the patient like a slide sheet and then inflated. The system is then inflated and the patient can be safely and effortlessly laterally transferred by two staff using the straps attached and involving no lifting. Easy to clean for infection control purposes. Can be used on the floor and with x-ray machines.

Porta 200



Max User Weight 30st

This portering Chair is excellent for hospital use, it features car type 3 point safety belt, independent pull-out leg rests, Drop side arms, Drip stand, sliding footplate and document pouch.

Seat width – 21”

Seat depth – 18” Seat height 21”

Bathroom Aids

Titan Shower Commode



Max User weight 30st

This Shower/commode chair features:

Flip away footplate

Flip away armrests

Padded Seat Cover

4 Braked castors

Optional drip stand, toilet roll holder, oxygen carrier

Seat Width – 46cm

Seat Depth – 46cm

Seat Height – 52cm

Titan Plus Shower Commode



Max User weight 50st

This Shower/commode chair features:

Flip away footplate

Flip away armrests

Padded Seat Cover – Bucket with lid

4 Braked castors

Optional drip stand, toilet roll holder, oxygen carrier

Seat Width – 54cm

Seat Depth – 43cm

Seat Height – 50cm

Zenith Bariatric Static Commode



Max User Weight 60st

This is a sturdy, static metal commode with potty and lid.

Features:

Detachable arms

Height adjustable

Seat Width – 20”

Seat Depth – 14”

Seat height – 18-24”

Rental Terms and Conditions

- 1. Delivery of a rental product can be secured by the health service with a signed rental agreement and purchase order number. Prices stated are ex vat and including delivery and re-collection.**
- 2. Private companies or individuals can secure rental by paying minimum rental period in advance.**
- 3. Emergency short term rentals available on request – delivery surcharge may be applicable.**
- 4. The minimum rental period for a product will be 4 weeks from date of delivery at a rate shown. These prices are exclusive of V.A.T. If the product is ready for return before the minimum rental period has expired you will still be charged the minimum rental. If rental exceeds 4 weeks payment due is per week or part of week thereafter.**
- 5. V.A.T. is chargeable at the standard rate.**
- 6. The chair will be delivered to your premises and the user instructed in its operation.**
- 7. We will require you to sign for the chair and to confirm that it was received in good condition. It is your responsibility to note any faults on receipt of product.**
- 8. When you have finished with the product you should fax us the cancellation form. RC-1 or telephone 028 92677077 to arrange collection. Rental Charges will only cease upon receipt of a signed cancellation form. We will then arrange collection of the product, where upon it will be checked for any damage**
- 9. The customer will be liable for any damage caused to the product, including upholstery damage and charged for any repair or replacement as appropriate.**

Purchase &/ Rental Prices available on request

John Preston & Co

7a Altona Rd

Lisburn

Antrim

BT27 5 QB

Key Contact - Lynne Breen Sales Manager

Direct line: 028 92633 790

lynne@johnpreston.co.uk

www.johnpreston.co.uk