



The ONLY bedpan system PROVEN to prevent the transmission of HCAs

So you can be sure *Quattro* is giving you the most advanced protection against the spread of infection, we partnered with TECcare® to develop a system that is proven to work, as verified by independent laboratories.

The safest and most effective macerator infection control system on the market.

Until now macerator suppliers have just sold disinfectant or added it to their machines' cycle to clean the hopper. We've gone one step further. Now every *Quattro* bedpan macerator has this pioneering system, able to deliver TECcare's outstanding antibacterial control. So not only can you depend on *Quattro* to run reliably. You can depend on it to protect against the spread of HCAs.

Introducing Haigh's Quattro, now with TECcare® CONTROL

TECcare® CONTROL

- ✓ Kills microbes within minutes of contact
- ✓ Effective on all microbial classes*
 - :: FRAGRANCE-FREE
 - :: CHLORINE-FREE
 - :: ALCOHOL-FREE
 - :: NON-CORROSIVE
 - :: NON-IRRITANT
 - :: NON-TOXIC

More than 50 NHS Trusts and a number of large care home groups now use TECcare® products – and this number is growing every day.

*MRSA and C. Diff account for less than 10% of HCAs so it's vital you use a broad-spectrum disinfectant like TECcare® CONTROL



Superior protection on the inside...

Quattro automatically injects TECcare® CONTROL into the lid spray to disinfect the hopper for up to 24 hours.

...and on the outside

Biomaster's silver antibacterial technology is built into *Quattro's* surfaces where it inhibits the growth of up to 99.99% of bacteria. This protection lasts throughout the macerator's working life.

Call us on **01989 760 200** or email sales@haigh.co.uk for more information.
We'll be happy to help you.



DESIGNED, MANUFACTURED AND ASSEMBLED IN THE UK BY HAIGH
RELIABILITY AS STANDARD

To achieve outstanding infection control you need to combine a number of things.

A reliable bedpan macerator must be one of them.

Human waste can easily be the cause of an HCAI as it can carry both viral and bacterial diseases. And as not all patients who carry a disease show signs of illness, every patient must be considered a potential risk. So the safe disposal of human waste must be at the heart of your infection control strategy.

In practice this means disposing of it quickly, with no bedpans piling up in your sluice room. And safely, using clean and reliable equipment. Reliable because a broken-down macerator not only fails to do its job, it presents a significant infection risk.

And you want, as far as possible, to eliminate that risk.

To do this you need *Quattro*.


4 outstanding levels of infection control

1 *Quattro* is the most reliable macerator on the market. We put it through thousands of repetitive operations to make sure it would stand up to the demands of the busiest wards. Which it does. Time after time after time.

2 NEW Only *Quattro* has the pioneering system able to deliver TECcare's outstanding antibacterial control.

3 *Quattro* has antibacterial protection from Biomaster built into its external surfaces.

4 We designed *Quattro's* smooth surfaces so they are easy to clean with no dirt traps where bacteria can hide and multiply.



:: *Quattro's* mechanical foot pedal is not dependant on an electronic sensor. So it will ALWAYS work.



:: You can use *Quattro's* soft touch button or the foot pedal to open the lid.

Call us on **01989 760 200** or email sales@haigh.co.uk for more information.

We'll be happy to help you.

