



A new benchmark for macerator infection control

We've collaborated with experts TECcare to create the only macerator with antimicrobial technology that is clinically proven. **Independent laboratory tests** confirmed Quattro used with TECcare CONTROL® gives 24 hours protection against all known hospital-based pathogens.

Haigh chose TECcare CONTROL as it's the **only** broad spectrum air and surface disinfectant effective against a wide range of pathogens that **doesn't corrode or damage the internal parts of the macerator.**

TECcare CONTROL is effective against all known pathogens including:

- | | |
|---|--|
| Acinetobacter Baumannii | Klebsiella Pneumoniae |
| Candida Albicans | MRSA Staphylococcus Aureus |
| Clostridium Difficile | Norovirus (Feline Calicivirus) |
| ESBL Escherichia Coli | Tetracycline Resistant Pseudomonas Aeruginosa |
| Hepatitis B+C Virus | Transmissible Gastroenteritis Virus |
| Human Immunodeficiency Virus type 1 (HIV 1) | Vancomycin Resistant Enterococcus Faecalis (VRE) |

..... CONTACT US FOR THE FULL LIST

It is specifically designed to reduce the risk of infection and cross contamination by interrupting key transmission pathways. By **cleaning and disinfecting simultaneously**, just a single stage process will achieve the highest level of cleanliness.

- | | | |
|---------------|----------------|-----------|
| Alcohol free | Fragrance free | Non-toxic |
| Chlorine free | Non-corrosive | |

Combined with the most reliable macerator, it provides you with the most effective and cost efficient infection control tool for disposing of pulp products.

Discover more at www.teccare.com or contact us:

T: +44 (0)1989 760 200

E: sales@haigh.co.uk



TECCARE® TECcare is a registered trademark of Talley Environmental Care Ltd