

Accessories

Wrist Worn Fall Detector (ZNC440)

- Pressure sensing technology used to detect auto cancel i.e. detects if the client has got back to their feet after a fall condition
- Solid state, dual axis accelerometer sensor
- Multi-directional movement and impact detection
- 4 hourly battery self-test and low battery reporting
- 5 sensitivity settings
- Reassurance vibration on fall detection
- 2-4 years battery life (typical - this will vary depending on number of calls)
- Alarm and event logging for reactive and preventative applications
- Can be worn on the wrist or on a belt




The wrist worn fall detector / pendant provides both a reactive and preventative solution that detects when the wearer has fallen, tripped or stumbled. The battery life is 2 - 4 years and will also report when the battery is getting low.

The fall detector uses both pressure sensing technology and a solid state, dual axis accelerometer to detect when a fall has occurred. It allows for further configuration by giving you the choice of 5 sensitivity settings so as to tailor the device to the individual. It is also equipped with 3 auto cancel modes as well as the ability to disable this facility.

Audio Visual Indicator (ZNC340)

- Compatible with CP3400 and CP5000 Call Processors
- Connected on a 2 pair power & data network
- 5 x high power Superflux RGB LED's
- 4 x high power Superflux white LED's
- Integral sounder: 80dB @ 1 metre
- Dimensions: 79 x 95 x 92mm (LxWxD)
- Weight: 164 grams



Our AVI3000 Audio Visual Indicator has been designed to meet the HTM 08-03 Colour Coding and Sound Pattern guidelines which makes it a perfect solution for an NHS Ward. The new model offers improved aesthetics, better light distribution making it easier to see, and an enhanced internal sounder. It is used above a door to highlight the location (room or bed), the type of call and where the call has been made from on the system and/or can be used as 'follow me' lights (commonly used in hospitals).

Bed Epilepsy Monitor (E40134)

- Detects nocturnal tonic-clonic seizures
- Adjustable duration/sensitivity
- No weight restrictions – suitable for adults and children
- Supplied with bed sensor mat
- Supplied with mains adapter and 2 x AA batteries
- Compatible with ATX3000plus, ATX4000/5000 and ATX3000 with adapter W07221



The Bed Epilepsy Monitor is designed for monitoring and detecting seizures with muscle jerking of a sleeping person with epilepsy. It can also give notification when a person leaves the bed and there are no more micro-movements, such as those caused by a person's breathing or heart beating. This can be especially useful for those who have a tendency to wander.

The Bed Epilepsy Monitor consists of two main components; a flexible and durable bed sensor placed under the mattress and a monitor with sophisticated embedded software. Together they detect the micro-movements of a person lying on the bed - even the heart beat and breathing, and the faster movements such as muscle spasms when person has a tonic-clonic seizure. When a person has faster and longer continued movements than pre-set thresholds the notification is triggered. This will sound an alarm (with adjustable volume) and send an "emergency" call via the ATX unit.

The kit includes a bed sensor, a wall mounted control unit, a plug-top power supply and a RJ11 lead for connection to an ATX3000 plus, ATX4000 or ATX5000 room unit. It can also be connected to older ATX3000 room units using the W07221 adapter lead.

Floor Pressure Mat (W07262)

- Dimensions - 92.5cm x 63.5cm

The Floor Pressure Mat is designed to combat the problem of residents who are inclined to walk undetected especially at night and are suffering falls and accidents.

It can be placed typically at the bedside or at the door threshold and plugged in to the call socket of your Nurse Call. It can then alert staff even before the resident is up and about, sounding the alarm with a small amount of pressure thus enabling staff to Investigate.

The Floor Pressure Mat has been proven to reduce the rate of falls and accidents considerably in numbers of Hospitals and Nursing Homes throughout the country. The Floor Pressure Mat has a heavy non slip backing, is professionally sealed so can easily be cleaned for liquid spills and is fully serviceable. The kit includes a RJ11 lead for connection to an ATX3000 plus, ATX4000 or ATX5000 room unit. It can also be connected to older ATX3000 room units using the W07221 adapter lead.



Tel: 0800 052 3616

Email: sales@aidcall.co.uk

Web: www.aidcall.co.uk  [@AidCall](https://twitter.com/AidCall)

Aid Call • Unit 4 • Brunel Court • Brunel Road • Newton Abbot • TQ12 4PB

"Aid Call" is a trading name of Tynetec Limited, one of the UK's leading suppliers of Telecare and remote health equipment.



Air Activated Switch (W07223)

- Flexible PVC actuator pad with 2 metre long PVC tubing.
- Free standing or fixed (mounting plate supplied)
- Actuator pad dimensions: 93mm diameter x 16mm deep



The air activated switch can be used in wet areas such as bath and shower rooms. The actuator pad is connected via PVC tubing to a wall mounted switch plate. The kit includes a RJ11 lead for connection to an ATX3000 plus, ATX4000 or ATX5000 room unit. It can also be connected to older ATX3000 room units using the W07221 adapter lead.

Bed & Chair Pressure Sensor (W07112 & W07114)

- Bed Sensor: 76cm x 28cm
- Chair Sensor: 38cm x 25.4cm

These are in-bed and in-chair devices that monitor occupancy, and will automatically raise alerts to carers when a potential risk situation is identified. This can be linked to light switching units that automatically switch on room lights when a client/patient leaves their bed during the night, providing dementia care facilities.

The kit includes a RJ11 lead for connection to an ATX3000 plus, ATX4000 or ATX5000 room unit. It can also be connected to older ATX3000 room units using the W07221 adapter lead.



PIR Sensor Kit & ATX5000 (ZNC293)

- Includes PIR kit with on/off switch & ATX5000 Call Point
- Switchable walk test LED
- Dual purpose optics: wide-angle 11x11m or long-range 1.7x17m in multi-level or pet alley patterns
- Selectable pulse (2 or 4)
- 2 minute interval period
- High detection sensitivity
- Tamper switch
- 10 year battery life (9V lithium battery - 1200mAH)
- Dimensions - H 114cm x W 70cm x D 76cm



Battery operated PIR Movement Detector for monitoring resident/patient activity or intruder detection. The kit includes an on/off switch as well as an ATX5000 Call Point.

Tel: 0800 052 3616

Email: sales@aidcall.co.uk

Web: www.aidcall.co.uk @AidCall

Aid Call • Unit 4 • Brunel Court • Brunel Road • Newton Abbot • TQ12 4PB
"Aid Call" is a trading name of Tynetec Limited, one of the UK's leading suppliers of Telecare and remote health equipment.

Big Red Button (W01060)

- Large momentary action push button
- Easily adjusted actuation force
- 15 settings from 200g to 1500g
- Dimensions: 130mm diameter x 40mm high
- Weight: 250 grams
- Lead length: 2 metres
- Additional adapter lead required (see options below)



The big button switch is a free-standing alarm button which can be placed on a table within reach of the user or screwed to the wall. The big button switch is suitable for persons with mobility or dexterity difficulties or for use as an easy activation panic button when connected to an ATX call point.

The force needed to activate the switch can be increased or decreased to suit the needs of the user. One press of the switch will initiate a call on the nurse call system.

Adapter Lead Options

The big button switch simply plugs into an ATX3000 plus, ATX4000 or ATX5000 nurse call room unit using one of the 3 adapter lead options listed below:

Part No.	Description
ZNC246	<ul style="list-style-type: none">• CALL adapter lead (2 metres)
ZNC254	<ul style="list-style-type: none">• EMERGENCY adapter lead (2 metres)
ZNC258	<ul style="list-style-type: none">• ASSISTANCE adapter lead (2 metres)

Tel: 0800 052 3616

Email: sales@aidcall.co.uk

Web: www.aidcall.co.uk  [@AidCall](https://twitter.com/AidCall)

Aid Call • Unit 4 • Brunel Court • Brunel Road • Newton Abbot • TQ12 4PB

"Aid Call" is a trading name of Tynetec Limited, one of the UK's leading suppliers of Telecare and remote health equipment.