

## The 'Nurse Alert Mat'

The 'Nurse Alert Mat' designed to combat the problem of residents who are inclined to walk undetected especially at night and are suffering falls and accidents.

The Mat can be placed typically at the bedside or at the door threshold and plugged into the call socket of the Nurse Call system or the mobile Floor Sentry Monitor. It can then alert staff even before the resident is up and about, sounding the alarm with a small amount of pressure thus enabling staff to Investigate.

It has also proven to reduce the rate of falls and accidents considerably in a number of Hospitals and Nursing Homes throughout the country and was featured in the Highland Regional Health Board Seminar as an extremely effective resident care product.

These mats can also be easily cleaned for liquid spills etc, is also fully serviceable and carries a 12 Month guarantee. The 'Nurse Alert Mat' is fitted with a 3mtr lead and comes with a stereo or mono jack plug for connection to a Standard Nurse Call socket. However, for a small additional cost we can supply or fit any length of lead, type of plug or adapter to suit your nurse call system.

For a Free No Obligation demonstration at your place of work then call Lotus Care on 0800 8499 121 and a member of our sales team will be glad to arrange a time convenient to you.



0 8 0 0 8 4 9 9 1 - 2 - 1

[www.lctuk.com](http://www.lctuk.com)

### The 'Nurse Alert Mat'

- 3M Cable
- Non-Slip High Grade Rubber
- Carpeted Waterproof Top
- Stereo/Mono Connector
- High Sensitivity Pressure Pad
- Sealed Waterproof Glue
- Correct Connector Supplied

